THE LOCOMOTOR RANGE™

Rota Cushion™ - Instructions

REF - MPC

LOCO-068	**Standard – 380mm diameter**
LOCO-191	**Large – 450mm diameter**



⚠ **All patient moving and handling tasks should be carried out in accordance with approved patient moving and handling techniques and after a full assessment of patient, equipment, task and environment. Always follow the approved policy of your place of work when adopting moving and handling practices**.

Description

The Rota Cushion is a soft padded, flexible turning disc which allows a smooth seated swivel transfer. It is ideal for use on car seats, chairs and beds.

Instructions for use

Place the Rota Cushion on car seat, chair, bed or similar. Sit centrally on to the Rota Cushion. Utilising the Rota Cushion rotate the whole body to the desired location. The swivel action of the Rota Cushion allows you to rotate smoothly into the correct position on the seated surface with minimum effort.

⚠ **For further information and training on moving and handling techniques please contact your Moving and Handling Advisor or similar**.

Care and use

⚠ Before each use inspect item. **Do not use** if the product is soiled or damaged

⚠ For seat transfers only. Do not use to stand on

⚠ Remove steel rotating disc before washing. Keep steel disc safe and replace when outer layers have been washed and dried. Do not submerge steel disc in water. If necessary, wipe steel disc clean with alcohol and dry thoroughly.

71° Washing temperature up to 71°C

Do not bleach

 Do not dry clean

 May be tumble dried on cool setting

 Maybe drip dried

Storage

 Store away from direct sunlight and sources of heat

 Store in dry conditions

---

Specifications

Material: 65% Polyester 35% cotton surface. 100% nylon lining. 5oz HSR poly T/L wadding.

Maximum user weight (MUW): 150 KG

 Not for lifting

 Does not contain latex

---

Disposal

 May be recycled where facilities exist

---

Warranty

 Warranted to specification for 12 months from date of purchase.

Applied regulations and standards

MD Class 1 Medical Device UK MDR 2002 EU MDR 2017/745 BS EN 12182: 2012 Tech Aids for disabled persons	CE UKCA

For further assistance please contact your supplier or customer services at:

Manufacturer

Select Healthcare (UK) Ltd.
Unit 2 Albert Road
Rushden
Northants
NN10 0BU

Tel: +44 (0) 1933 411729
Fax +44 (0) 1933 410213
Email: info@selecthealthcare.co.uk
Website:
www.selecthealthcare.co.uk

European Representative

EC REP **Select Healthcare (EU),** 72 Liffock Park, Coleraine, Northern Ireland, BT51 4DG

Glossary of symbols and terms used in product and packaging labels that are not expressed above.

Symbol	Meaning	Associated Terms
	Consult instructions	
DI	Device Identifier – (GS1 UK GTIN 13 used)	AI (01) – Application identifier 01
SN	Serial Number	AI (21) – Application identifier 21
LOT	Lot / Batch Number	AI (10) – Application Identifier 10
	Date of manufacture	AI (11) – Application Identifier 11
UDI	Unique Device Identifier - collective term for DI, SN, LOT & date of manufacture (all contained within machine readable GS1 DataMatrix)	GS1 DataMatrix = 2D Matrix symbols that can carry all GS1 keys and attributes